






CLOSED MODULAR TROLLEY MMSlog®

MAIN CHARACTERISTICS OF THE CLOSED MODULAR TROLLEY

- Closed Modular Trolley for transportation and storage of healthcare material within the medical center, dimensions in accordance with ISO 3394.
- Remarkable features: modularity (available in two different heights and with one or two body structures), lightness, stability, robust construction, easy to handle, clean and sanitize, perfect for the surgical area.
- Trolley body cast in one piece of folded steel-sheet with epoxy coating (in color blue RAL 5024*).
- The trolley upper surface is a 10 mm thick PVC in solid white color.
- Trolley base with bumpers (ABS) to absorb shocks, 4 double-rolling band rotary wheels (diam. 125 mm), two of them with brake. Optionally wheels can be with double band.
- Front side closing with roller-doors (ABS), made in one piece and removable, easy opening and with a lock in the upper part.
- Both sides of the trolley are set up with two ergonomic handles, for an easier manipulation.
- Inside body equipped with sliding-rails side panels designed for 600*400 mm modules (ISO3394) and also for other storing accessories.

REFERENCE INDEX AND CHARACTERISTICS OF THE AVAILABLE MODELS

Image	Reference	Description	External Dimensions mm.			Weight
			Width	Depth	Height	Kg.
	TC101060151	Closed Modular Trolley – 1 body front 600 mm. And height 1500 mm.	830	570	1670	65
	TC101060181	Closed Modular Trolley – 1 body front 600 mm. and height 1800 mm.	830	570	1960	70
	TC201060151	Closed Modular Trolley – 2 bodies front 600 mm. and height 1500 mm.	1545	570	1670	125
	TC201060181	Closed Modular Trolley – 2 bodies front 600 mm. and height 1800 mm.	1545	570	1960	135
	TC101040151	Closed Modular Trolley – 1 body front 400 mm. and height 1500 mm.	630	770	1670	70
	TC201040151	Closed Modular Trolley – 2 bodies front 400 mm. and height 1500 mm.	1145	770	1670	128
	SW594012504	Set of 4 wheels (2 with brakes) of double band for Modular Trolley MMSlog®	Ø 125	-	150	-

(*) The standard colour available is RAL 5024. Nevertheless, for orders equivalent to 25 or more body trolleys, they can be manufactured in other RAL colours.