



MMSlog® OPENED TRANSPORT TROLLEY (W/ MIX MODULES & SHELVES)

MMSlog® OPENED TRANSPORT TROLLEY, 1 X 10 GUIDE RAILS + 1 X 3 SHELVES (CT431000000)

- Opened trolley for transport & picking of sanitary products within the health center
- Main features: easy to handle, light-weight, stable, easy to clean & sanitize
- Body cast in one piece of steel tube 40*20 mm with epoxy coating
- Equipped with 1 column of 10 pairs of ABS sliding guide rails with “stop” reed at the ends
- The 1st column is set for the use of 600*400 mm standard modules**
- The 2nd column is fitted with 3 HPL shelves of dimensions 850/420*605*10 mm to carry boxes
- The 4 edges are provided with 1 ergonomic handle for easy handling
- Rubber bumpers on the base
- 4 anti-thread rotary wheels, 2 with brake, diameter 125 mm

Reference	Size			Núm.	Pairs	Casters
Code	Height	Width	Depth	Bodies	Guides/shelves	Diameter – Brakes
CT431000000	1555 mm	1450 mm	745 mm	1+1	10+3 shelves	125 mm - 2 c/brake